

OVERVIEW

The H12775 is a photon counting head device containing of a 10-mm (1/2") diameter head-on photomultiplier tube, along with a high-speed photon counting circuit and a high-voltage power supply circuit. The high voltage supply for photomultiplier tube and the discrimination level are preset to optimum values, allowing photon counting measurement by just connecting a +5 V supply. The H12775 has a function of over light detection output for indicating the excessive light operating condition.



PRODUCT VARIATIONS

Type No.	Spectral response	Photocathode	Window material	Features
H12775	300 nm to 650 nm	Bialkali	Borosilicate glass	For visible range

This product can't be used at vacuum environment or reduced pressure environment.

SPECIFICATIONS

(at +25 °C)

Parameter		Description / Value		Unit
Input voltage		+4.75 to +5.25		V
Max. input voltage		+6		V
Max. input current		40		mA
Effective area		φ10		mm
Peak sensitivity wavelength		420		nm
Count sensitivity	300 nm	Typ.	2.5×10^5	s ⁻¹ ·pW ⁻¹
	400 nm	Typ.	4.4×10^5	
	500 nm	Typ.	3.6×10^5	
	600 nm	Typ.	1.2×10^5	
Count linearity *1		5.0×10^6		s ⁻¹
Dark count *2		Typ.	30	s ⁻¹
		Max.	100	
Pulse-pair resolution		20		ns
Output pulse width		10		ns
Output pulse height	Load resistance 50 Ω	Min.	+2.0	V
		Typ.	+2.2	
	Un-terminated	Min.	+4.0	
		Typ.	+4.4	
Recommended load resistance		50		Ω
Signal output logic		Positive logic		—
Over light detection output *3	High level	Min.	+3.5	V
	Low level	Max.	+0.5	
Operating ambient temperature *4		+5 to +40		°C
Storage temperature *4		-20 to +50		°C
Weight		79		g

*1: Random pulse, at 10 % count loss *2: After 30 minutes storage in darkness *3: Load resistance 10 kΩ *4: No condensation

PHOTON COUNTING HEADS H12775

Figure 1: Count sensitivity




Figure 2: Dark count




Figure 3: Output waveform




Figure 4: Count rate linearity and over light detection output





Figure 5: Count rate linearity correction



Linearity correction formula


$$N = \frac{M}{1 - Mt}$$

N: real count rate (s⁻¹)

M: measured count rate (s⁻¹)


t: pulse pair-resolution (s)

Figure 6: Block diagram



TPMOC0264EC

Figure 7: Dimensional outline (Unit: mm)



* Option: Available with BNC/SMA connector.

TPMOA0101EA

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: info@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 inf. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jianning Center, 27 Dongshanhuan Beilu, Chaoyang District, 100020 Beijing, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: Hamamatsu Photonics Taiwan Co., Ltd.: 8F-3, No. 158, Section 2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)03-659-0080, Fax: (886)03-659-0081 E-mail: info@tw.hpk.co.jp

TPMO1084E01

OCT. 2017 IP